# SOL'S WILSON BW WOMEN



## WOMEN'S LIGHTWEIGHT BODYWARMER

02890











**POLYAMIDE 380T** Zip opening 2 front zip pockets 2 inside pockets Small high collar Close-fitting cut "In bag" folding system with stoppers

#### Care





## Composition

Outer fabric and lining: 100% polyamide. Padding (please note that we are in the process of transitioning from stock): 90% down and 10% feather until stock lasts. Then transition to 60% Sorona 40% recycled polyester.

### **Available colours**











## A committed product



## **Related products**





SOL'S WILSON WOMEN 02899


## SOL'S WILSON BW WOMEN



WOMEN'S LIGHTWEIGHT BODYWARMER

02890

#### **Available sizes**



## **Packaging**

Carton size :  $60 \times 40 \times 25 \text{ cm}$ 

Weight per carton: 3.00 kg

自 10 📙 1

#### **Customisation**

- Embroidery: This technique is generally used for personalisations that aim for a high-end finish. This technique is the most resistant to washing and use. The embroidery can be applied directly to the product or through embroidered patches. It can be done with thickness effects (with foam) or through patches that will then be affixed to the final product, allowing for variations in materials.
- Flex: It is the recommended marking technique for small and medium-sized series. There are several aspects: smooth, velvet, fluorescent, glitter, gold and silver. These are iron-on vinyls that are cut and glued by hot pressing. They are very well suited to a wide range of materials and substrates.
- Transfert: The right technique for all materials. It is recommended for luggage, heavy garments, difficult to access surfaces. It consists of transferring the marking from a paper support to the garment by hot gluing. The marking by gluing brings rigidity to light supports at the level of the marking zone, nevertheless, the product keeps all its qualities of comfort.